ARTHROSAMID<sup>®</sup>

# Arthrosamid®

<b>en</b> .....	<b>4</b>
<b>de</b> .....	<b>7</b>
<b>fr</b> .....	<b>10</b>
<b>es</b> .....	<b>14</b>
<b>da</b> .....	<b>17</b>
<b>fi</b> .....	<b>20</b>
<b>it</b> .....	<b>24</b>
<b>nl</b> .....	<b>27</b>
<b>no</b> .....	<b>31</b>
<b>pl</b> .....	<b>34</b>
<b>pt</b> .....	<b>37</b>
<b>sv</b> .....	<b>41</b>
<b>tr</b> .....	<b>44</b>

## INSTRUCTIONS FOR USE

### PRODUCT DESCRIPTION

Arthrosamid® is a non-absorbable, biocompatible, injectable, transparent hydrogel for intra-articular injection. Arthrosamid® consists of 2.5 % cross-linked polyacrylamide and 97.5 % non-pyrogenic water. The hydrogel is supplied in a sterile, prefilled 1 ml syringe sealed with a Luer Lock fitting and a tip cap.

Arthrosamid® is intended to be injected intra-articularly with a sterile 21Gx 2 inches (0.8x50mm) needle with Luer Lock fitting. A 3-part label with the LOT number and expiry date is available on the blister pack. Attach one label to the patient record to ensure product traceability.

### INDICATIONS, PATIENT GROUP AND USAGE

Arthrosamid® is intended to be used for symptomatic treatment of adult patients with knee osteoarthritis.

### MODE OF ACTION

- Arthrosamid® integrates into the synovial tissue of the inner joint capsule and decreases joint stiffness
- Arthrosamid® forms a cushion-like membrane comprising the hydrogel incorporated in a matrix of fine tissue network covered by a synovial lining
- Arthrosamid® diminishes pain and improves function of the knee affected by osteoarthritis

The non-absorbable, non-bio-degradable and non-migratory characteristics of Arthrosamid® provide cushioning of the inner joint capsular tissue.

### CONTRAINDICATIONS

Arthrosamid® should not be injected:

- If an active skin disease or infection is present at or near the injection site
- If the joint is infected or severely inflamed
- If a degradable intra-articular injectable such as hyaluronic acid is present, it must be expected to be absorbed according to manufacturer's information for the specific product, before injection with Arthrosamid®
- If the patient has previously received treatment with a different non-absorbable injectable/implant
- If the patient has received a knee alloplasty or has any foreign material in the knee
- If the patient has undergone knee arthroscopy within the last 6 months
- In hemophilia patients or in patients in uncontrolled anti-coagulant treatment

Do not subsequently inject a different non-absorbable implant.

### WARNINGS

- Do not inject intravascularly, extra-articularly or in the synovial tissue of the capsule
- Do not inject corticosteroids along with Arthrosamid®

### PRECAUTIONS

Usual precautions associated with invasive joint procedures should be followed.

Patients with acute or chronic infections in other sites of the body should be treated with caution.

Arthrosamid® should be used with caution in patients with e.g. autoimmune disorders and in patients with uncontrolled diabetes as well as in patients undergoing major dental work or surgery.

Safety and effectiveness has not been established in patients under 18 years, in pregnant or lactating women, or in patients having foreign material implanted in the knee.

Only use the product if the packaging and products are intact and undamaged. Do not re-sterilize Arthrosamid®.

Arthrosamid® is only intended for use as an intra-articular injectable. As with any invasive joint procedure, there is a small risk of infection when injecting Arthrosamid® and therefore a single dose of prophylactic antibiotics must be administered prior to injection.

The Arthrosamid® syringe is intended for single patient use – do not store unsealed syringes and re-use. Re-use increases the risk of contamination and thereby increases the risk of infection. Do not use Arthrosamid® once expired.

Do not inject any pharmaceuticals or biological substances into the hydrogel.

### **PROPHYLACTIC ANTIBIOTIC TREATMENT**

Arthrosamid® must be administered with prophylactic antibiotics. The following combination of antibiotics is recommended: Azithromycin 500 mg p.o. and Moxifloxacin 400 mg p.o. administered 1-6 hours prior to injection.

The above combination of antibiotics will reach a high concentration in the tissue at the time of injection and should be given only once. This combination covers up to 95 % of the normal occurring skin flora (both aerobic and anaerobic species) and has a long half-life.

### **METHOD OF ADMINISTRATION**

Arthrosamid® must be administered by qualified physicians familiar with intra-articular injection procedures, such as orthopaedic surgeons or rheumatologists. Inject into the joint space only.

The procedure must be conducted under aseptic conditions and using no-touch technique. It is essential that at least 5 cm around the injection site is swabbed prior to injection using e.g. chlorhexidine with alcohol three times with one minute interval. The injection of the hydrogel must be performed under local anaesthesia.

Remove joint effusion, if present, before injecting Arthrosamid®, using the same needle. The recommended needle is standard 21G x 2 Inches (0.8x50 mm) with Luer Lock.

Recommended dosage: 6 ml of Arthrosamid®.

Injection should be performed with the help of ultrasound in the lateral, proximal recess. If ultrasound is not available, the injection can be performed along the joint line lateral to the ligamentum patellae. The patient can be in a supine position or sitting with the knee bent between 80°- 90°.

Remove the protective Tip Cap from the Arthrosamid® syringe. Secure the needle firmly into the Luer Lock socket and make sure the needle is mounted correctly. In the event of an obstruction of the needle, do not increase the pressure on the plunger. Stop the injection and replace the needle.

The syringe, the needle and any unused material must be discarded after the treatment session. The needles are discarded in a disposal container according to normal practice/local regulation. Used syringes may be discarded as normal waste.

### **PATIENT INFORMATION**

Please inform your patient

- As with any invasive joint procedure it is recommended to avoid strenuous activity (e.g. tennis, jogging or long walks) during the first few days after the injection.
- Advise patients that if they require a surgical or dental procedure in the future, they should tell the treating physician that they have a permanent implant and they should discuss the possible need for prophylactic antibiotic therapy with their treating physician.

Patients should inform health care professionals of their Arthrosamid® implant for accurate future medical assessment.

The patient should be informed about indications, expected results, contraindications, warnings, precautions and potential complications.

In case of complications, the patient should contact the injecting physician immediately for treatment.

## ADVERSE EVENTS

Common adverse events:

According to a clinical study, injection-related mild to moderate pain and/or mild swelling during the first weeks to months following the injection may occur. Swelling of Baker's cyst may also occur.

Rare adverse events:

Infections might occur at the injection site and must be treated immediately with antibiotics.

In case of infection, the use of corticosteroids should be avoided as this may worsen the bacterial infection and complicate its treatment.

## REPORTING OF ADVERSE EVENTS

All adverse events/complications must be reported to the local distributor or directly to Contura International A/S via Email to: [complaints@contura.com](mailto:complaints@contura.com).

## SUPPLY AND STORAGE

Sterile: Arthrosamid® is supplied sterile. The hydrogel must be kept away from sunlight. Do not freeze. The expiry date is indicated on the package.

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